Opis

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Instructions for use of the procedure pack

Description

The procedure pack designed for separation and concentration of bone marrow aspirate to obtain concentrated bone marrow fraction (BMAC). BMAC is designed for local injections to improve the regenerative properties of cells and tissues by providing a large number of stem cells and growth factors to specific areas.

Procedure pack consists of the following medical devices:

- ① BMAC, *Biovico* (separation device), 3 pcs.
- ② Microlance™ 3, 21G x 1½", *BD* (needle), 1 pc.
- ③ Microlance™ 3, 18G x 1½", *BD* (needle), 1 pc.
- ④ dicoNEX Single use syringe (with Luer-lock), 3-part 10 ml, *Zarys*, 4 pcs.
- ⑤ KD-STOP, Combi stopper, *KD Medical* (syringe cap), 1 pc.

Cautions

To be used exclusively by qualified, trained personnel.
Before using the procedure pack, please read the instructions for use.
The procedure pack consists of disposable medical devices. Do not reuse.
Do not use any of the medical devices in case of damage of its direct package.
Medical devices are sterile only during their shelf life indicated on the package.
Store the procedure pack at a temperature between +5°C and +25°C. Protect from direct sunlight, high temperature and high humidity. Keep out of the reach of children.
Wear personal protective equipment when using the procedure pack.
Ensure that no one is exposed to direct contact with collected biological material.
Place all used pack parts in a medical waste container and utilize them.
Local symptoms such as pain, redness and swelling may occur at the injection site.

Contraindications

Children, pregnant or nursing women should not be treated with the procedure pack as there are no clinical data available on therapy in those groups of patients.

NOTE: Use only freshly collected bone marrow for isolation. Isolated bone marrow concentrate should be used immediately after preparation. The bone marrow aspirate should be collected into the syringe filled with anticoagulant. In case of 3.13% sodium citrate, the volume ratio of anticoagulant to bone marrow aspirate should be 1:9. In order to obtain the best quality preparation, it is recommended that the bone marrow concentrate is prepared in accessories rinsed with anticoagulant. To avoid contamination of the marrow concentrate, perform the entire isolation procedure with caution to sterility.

Instructions for use

1. Remove all three BMAC tubes ① from the packages, position them vertically and mark with the patient's code or name to avoid confusion.
2. Check if the bottom screw cap of each BMAC tube ① is tightened correctly. Insert the needle ② into the center of the

- top cap and carefully transfer the marrow aspirate with anticoagulant into the BMAC tubes ①, filling each of them to 15 ml. The tubes should be placed vertically on a flat surface.
3. Open the centrifuge cover, place the BMAC tubes ① and properly prepared counterweights inside (if applicable).
4. Centrifuge for **15 minutes at 600 x g** (RCF).
5. Open the cover and carefully take out each BMAC tube ① with the centrifuged bone marrow aspirate. Hold BMAC tubes ① upright, do not shake them in order to avoid mixing of the separated fractions.
6. Replace the bottom screw cap of the first BMAC tube ① with the controller screw cap and then fasten the controller's bolt to the BMAC tube ①.
7. Turn the bottom part of the controller's bolt counterclockwise, until the plasma reaches the last ring in the top part of the BMAC tube ①.

NOTE: Each time before connecting the syringe to the BMAC tube ①, check if the syringe plunger is moving properly.

8. Open the top cap of the BMAC tube ①.
9. Attach the 10 ml syringe ④ in place of the top cap to collect platelet-poor plasma (PPP). Rotate the controller's bolt until the red blood cells (RBC) fraction reaches the line marked as '1' on the BMAC tube ①. Disconnect the syringe ④ with PPP. Then, to collect 2 ml of BMAC, connect another 10 ml syringe ④ in place of the top cap and turn the controller screw until the RBC fraction is at the level of the last ring in the top of the BMAC tube ①. Collect 1 ml of the large stem cell fraction by filling the syringe ④ up to the 2 ml, according to the scale on the syringe.
10. Disconnect the syringe ④ with the collected ready-to-use BMAC and close tightly with the syringe cap ⑤.
11. Mark the syringe ④ with collected BMAC with the patient's code or name to avoid confusion.
12. Repeat steps 6-10 for the second and third BMAC tubes, collecting BMAC into the same 10 ml syringe ④ used to collect the marrow concentrate from the first BMAC tube, subsequently filling it up to the 4 and 6 ml level.
13. After collecting the third fraction of BMAC into the 10 ml collection syringe ④, connect the needle ②. The bone marrow concentrate is ready for injection.

Organization assembling procedure pack:

Biovico Sp. z o.o.
Hutnicza 15 B
81-061 Gdynia, Poland, EU
+48 58 660 44 88
office@biovico.com
www.biovico.com

Date of issue: 07/2024 Rev. 01/2024

BIOVICO
MEDICAL BIOTECH