



SIRONIX | Healthium

# STATIV®



**ALL SUTURE ANCHOR | ROTATOR CUFF REPAIR**  
(2.5mm/2.9mm)



**EU-MDR**  
2017/745

**USFDA**  
510(K) CLEARED

**MAY THE SPORT OF LIFE NEVER STOP®**

# SCIENTIFIC EVIDENCE

## Research Article

A multicentre, retrospective and observational study to evaluate safety and functional outcomes of arthroscopic shoulder ligament repair using Sironix suture anchor.

**Authors:**  
Dr. K. Vijayaraj and Dr. V. Rathika.



## Research Article

Assessing clinical and patient reported outcomes of Sironix suture anchors in rotator cuff and Bankart repair surgeries.


**Authors :**  
Dr. Puthumanapully Vinod Kumar, Dr. Subin Sugath and Dr. Vijay Mohan.




# AVAILABLE OPTIONS

Self Tap	2.5mm	2.9mm
Suture	Double	Double / Triple
Tape	Double	Double
Drill Size	2.5 mm	2.9 mm

Note : Suture Size USP 2; Tape Size 1.5mm




**STATIV®** features a patented all-suture design that delivers secure fixation with minimal bone removal.



### Self-Tapping Anchor\*

Enables easy deployment while reducing the need for additional instrumentation.



### Tripod Bunching Design\*

Unique bunching pattern for secure cortical fixation.

### UHMWPE Suture Construct\*


High-strength, durable sutures for reliable soft tissue repair.

### Different Suture/ Tape Color option\*

Blue and black suture/tape for clear identification and easier suture management.

### Wide Range\*

Multiple suture/tape options to suit varied surgical techniques.



\*Data on file

# ORDERING INFORMATION

Code	Description
S33-2502-ST	Knotted UHMWPE Suture Self Tap Anchor-2.5mm Double loaded sutures-USP #2
S33-2503-ST	Knotted UHMWPE Suture Self Tap Anchor-2.5mm Triple loaded sutures-USP #2
S33-25E2-ST	Knotted UHMWPE Suture Self Tap Anchor-2.5mm Double loaded tapes-1.5mm
S33-2902-ST	Knotted UHMWPE Suture Self Tap Anchor-2.9mm Double loaded sutures-USP #2
S33-2903-ST	Knotted UHMWPE Suture Self Tap Anchor-2.9mm Triple loaded sutures-USP #2
S33-29E2-ST	Knotted UHMWPE Suture Self Tap Anchor-2.9mm Double loaded tapes-1.5mm

## ASSOCIATED INSTRUMENTS

Code	Description
S45-0029-LST	Stativ Sleeve For Self Tap 2.5mm/2.9mm
S45-ST25-DU	2.5mm Stativ Anchor Drill
S45-ST30-DU	2.9mm Stativ Anchor Drill

# SIRONIX

Scan the QR code to visit our website



STAY  
CONNECTED



Sironix Medtech



Sironix-Medtech



@sironix\_arthroscopysolutions



teamsironix@healthiummedtech.com

## MAY THE SPORT OF LIFE NEVER STOP®



Healthium Medtech Limited

#472-D, 13th Cross, 4th Phase, Peenya Industrial Area, Bangalore - 560 058. India

Tel: +91 80 4186 8198 | Email: care@healthiummedtech.com

Disclaimer:

Responsibility for administering products mentioned and interpretation of clinical data provided lies with the registered medical practitioner. For more information about products, please visit our website and refer to Instructions for Use document (IFU) document.

Reference: \*Data on file

SK005711082025|Stativ RC Brochure\_D|5000|DIGI